

Code No.19-0727

## SPECIFICATION

pH4.30 Sodium Citrate Buffer Solution (5 times concentrated)  
for Amino Acid Automated Analysis

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
* pH at 25°C	4.27~4.33
Concentration of Cl (Cl <sup>-</sup> )	0.87~1.07mol/l
Concentration of Na (Na <sup>+</sup> )	2.10~2.40mol/l
* Ninhydrin reactive substances	to pass test

\*:Specification of 5 times diluted solution