

Code No.19-0743

SPECIFICATION

Sodium linear -Dodecylbenzenesulfonate Standard for ABS Analysis

REQUIREMENT	SPECIFICATION
Appearance	White, powder or mass
Loss on drying at 105°C	max.0.5%
pH of a 10g/L solution at 25°C	5.0~7.0
Sodium sulfate (Na ₂ SO ₄)	max.0.1%
Dodecylbenzene (C ₁₂ H ₂₅ C ₆ H ₅)	max.0.1%
High performance liquid chromatography test	to pass test
Assay (mixture of isomers)	min.99.0%