

## SPECIFICATION

## Sodium Azide

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (NaN <sub>3</sub> ) Solubility in water Alkalinity (as NaOH) Loss on drying at 105°C Chromium (Cr)	White ~ nearly white, crystalline powder or mass min.98.0%(mass/mass) to pass test max.0.2%(mass/mass) max.0.1%(mass/mass) max.2ppm(mass/mass)

JIS K 9501 : 2019(\*:Additional test performed by Wako but not required by JIS)