

## SPECIFICATION

### 10% Silver Nitrate-impregnated Silica Gel

for Dioxins Analysis

REQUIREMENT	SPECIFICATION
Appearance	White small granules
Particle size (75~150 $\mu$ m)	min.75%
Particle size (pass 45 $\mu$ m)	max.5%
Precipitation volume	1.8~2.7ml/g
Loss on drying at 105°C	1.0~2.0%
Practical test	to pass test
Silver nitrate (AgNO <sub>3</sub> ) (calculated on the dried basis)	8~12%