

Code No.19-1397

SPECIFICATION

Sodium L-Tartrate

Japanese Pharmaceutical Excipients

REQUIREMENT	SPECIFICATION
* Description	Colorless crystals or white crystalline powder
Identification	to pass test
Optical rotation [α]D ₂₀	+25.0~+27.5°
pH	7.0~9.0
Purity	-
(1) Clarity and color of solution	to pass test
(2) Sulfate	not more than 0.019%
(3) Oxalate	to pass test
(4) Heavy metals	not more than 10ppm
(5) Arsenic	not more than 2ppm
Loss on drying	14.0~17.0%
Assay (as anhydrous) (after drying)	not less than 98.5%

Japanese Pharmaceutical Excipients(*:Reference item)