

## SPECIFICATION

## Salicylamide

for the Japanese Pharmacopoeia General Tests

REQUIREMENT	SPECIFICATION
Description (JP method)	White crystals or crystalline powder
Melting point (JP method)	139~143°C
Purity (JP method)	-
Ammonium (JP method)	to pass test
Loss on drying (JP method)	max.0.5%
Residue on ignition (JP method)	max.0.1%
Assay (JP method)	min.98.5%
Solubility in ethanol	to pass test
Assay (HPLC)	min.99.0%