

SPECIFICATION

Saikosaponin b2

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and Thin-Layer Chroma

REQUIREMENT	SPECIFICATION
Description (JP method)	White crystals or crystalline powder
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (JP method)	to pass test
Unity of peak (JP method)	to pass test
Assay (qNMR) (JP method)	under investigation
* Solubility in methanol	to pass test
* Loss on drying over silica gel in vacuum	max.2.0%

This product is suited for Saikosaponin b2 for assay (for assay 2) and thin-layer chromatography of general tests on The Japanese Pharmacopoeia. (*:Additional test performed by Wako)