

## SPECIFICATION

## Sodium Chloride

TraceSure(R)

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (NaCl)	min.99.95%(mass/mass)
Water-insoluble matter	max.0.005%(mass/mass)
pH of a 50g/l solution at 25°C	5.0~7.5
Bromide (Br)	max.0.05g/kg
Iodide (I)	max.0.05g/kg
Nitrate (NO <sub>3</sub> )	max.5mg/kg
Phosphate (PO <sub>4</sub> )	max.1 mg/kg
Sulfate (SO <sub>4</sub> )	max.0.01g/kg
Silicate (as Si)	max.0.02g/kg
Potassium (K)	max.0.03g/kg
Magnesium (Mg)	max.3mg/kg
Calcium (Ca)	max.0.01g/kg
Barium (Ba)	max.0.01g/kg
Copper (Cu)	max.1 mg/kg
Lead (Pb)	max.1 mg/kg
Iron (Fe)	max.0.5mg/kg
Ammonium (NH <sub>4</sub> )	max.5mg/kg

JIS K 8005:2016

\*:Additional test performed by Wako but not required by JIS