

# SPECIFICATION

## Sodium Hydroxide for Nitrogen Compounds Analysis for Nitrogen Compounds Analysis

REQUIREMENT	SPECIFICATION
* Appearance	White shot
Assay (NaOH)	min.96.0%(mass/mass)
Solubility in water	to pass test
Nitrogen compounds (as N)	max.5ppm(mass/mass)
Iron (Fe)	max.5ppm(mass/mass)

JIS K 8826:2020 (\*:Additional test performed by Wako but not required by JIS)