

Code No.20-0102

# SPECIFICATION

## Mercaptoacetic Acid

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless ~ nearly colorless, clear liquid
Assay (HSCH <sub>2</sub> COOH)	min.90.0%(mass/mass)
Solubility in water	to pass test
Residue after ignition (as sulfate)	max.0.01%(mass/mass)
Suitability for iron determination	to pass test
Iron (Fe)	max.5ppm(mass/mass)

JIS K 8630:2019(\*:Additional test performed by Wako but not required by JIS)