

Code No.20-0148

## SPECIFICATION

## Tin

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
Appearance	Silver gray shot
Copper (Cu)	max.0.005%(mass/mass)
Lead (Pb)	max.0.01%(mass/mass)
Antimony (Sb)	max.0.01%(mass/mass)
Bismuth (Bi)	max.0.003%(mass/mass)
Iron (Fe)	max.0.005%(mass/mass)

JIS K 8580:2011

\*:Additional test performed by Wako but not required by JIS  
This product is shot