

## SPECIFICATION

## Tin(II) Chloride Dihydrate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (SnCl <sub>2</sub> ·2H <sub>2</sub> O) Solubility in hydrochloric acid Sulfate (SO <sub>4</sub> ) Sodium (Na) Potassium (K) Copper (Cu) Calcium (Ca) Lead (Pb) Arsenic (As) Iron (Fe) Ammonium (NH <sub>4</sub> )	White crystals min.97.0%(mass/mass) to pass test max.0.002%(mass/mass) max.0.005%(mass/mass) max.0.002%(mass/mass) max.0.001%(mass/mass) max.0.002%(mass/mass) max.0.003%(mass/mass) max.2ppm(mass/mass) max.0.002%(mass/mass) max.0.002%(mass/mass)

JIS K 8136:2017\*:Additional test performed by Wako but not required by JIS