

Code No.20-0158

## SPECIFICATION

## Tin(II) Chloride

Wako 1st Grade

| REQUIREMENT                         | SPECIFICATION   |
|-------------------------------------|---|
| Appearance                          | White ~ nearly white,<br>crystalline powder or granules<br>to pass test |
| Solubility in dil.hydrochloric acid |   |
| Sulfate (SO <sub>4</sub> )          | max.0.03%   |
| Iron (Fe)                           | max.0.005%  |
| Assay                               | min.97.0%   |