

SPECIFICATION

Trifluoroacetic Acid for Amino Acid Sequence Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (50vol%, 260nm)	max.1.0
Absorbance (50vol%, 270nm)	max.0.10
Absorbance (50vol%, 280nm)	max.0.02
Absorbance (50vol%, 300~400nm)	max.0.01
Aldehydes (as CH ₃ CHO)	max.5ppm
Suitability for amino acids sequence analysis	to pass test
Assay	min.99.0%