

Code No.20-0706

SPECIFICATION

1,1,2,2-Tetrachloroethane

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless ~ nearly colorless, clear liquid
Assay (CHCl ₂ CHCl ₂)(GC)	min.97.0%(mass/mass)
Density (20°C)	1.593~1.601g/ml
Refractive index 20°C	1.493~1.497
Water	max.0.1%(mass/mass)
Acidity (as HCl)	max.0.003%(mass/mass)

JIS K 9520:2013

(*:Additional test performed by Wako but not required by JIS)