

SPECIFICATION

2,2,4-Trimethylpentane

for Beer Bitterness Analysis

| REQUIREMENT | SPECIFICATION |
|---------------------------|------------------------|
| Appearance | Colorless clear liquid |
| Density (20°C) | 0.690~0.693g/ml |
| Refractive index 20°C | 1.390~1.393 |
| Absorbance (275nm) | max.0.005 |
| Water | max.0.01% |
| Residue after evaporation | max.0.001% |
| Assay (capillary GC) | min.99.7% |