

SPECIFICATION

Tetrahydrofuran, Dehydrated, Stabilizer Free for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.884~0.889g/ml
Refractive index 20°C	1.406~1.409
Absorbance (250nm)	max.0.50
Absorbance (280nm)	max.0.20
Absorbance (350nm)	max.0.02
Water	max.0.005%
Peroxides (as H ₂ O ₂)	under investigation
Assay (capillary GC)	min.99.5%