

## SPECIFICATION

Tetrahydrofuran, Stabilizer Free  
for Spectrochemical Analysis

| REQUIREMENT            | SPECIFICATION          |
|------------------------|------------------------|
| Appearance             | Colorless clear liquid |
| Density (20°C)         | 0.885~0.888g/mL        |
| Refractive index 20°C  | 1.406~1.409            |
| Absorbance (215nm)     | max.0.50               |
| Absorbance (230nm)     | max.0.40               |
| Absorbance (240nm)     | max.0.30               |
| Absorbance (254nm)     | max.0.20               |
| Absorbance (270nm)     | max.0.05               |
| Absorbance (280nm)     | max.0.03               |
| Absorbance (300~400nm) | max.0.01               |
| Water                  | max.0.02%              |
| Fluorescence test      | to pass test           |
| Assay (capillary GC)   | min.99.9%              |