

## SPECIFICATION

0.1 vol% Trifluoroacetic Acid-Acetonitrile  
for HPLC

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (200nm)	max.1.00
Absorbance (210nm)	max.1.50
Absorbance (220nm)	max.1.50
Absorbance (230nm)	max.0.50
Absorbance (240nm)	max.0.30
Absorbance (254nm)	max.0.10
Absorbance (260nm)	max.0.03
Absorbance (270~400nm)	max.0.01
Water	max.0.05%
Fluorescence test	to pass test
Trifluoroacetic acid	0.18~0.20%