

Code No.20-1660

## SPECIFICATION

## Tetrahydrofuran

for RoHS

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Assay (C <sub>4</sub> H <sub>8</sub> O) (GC)	min.99.5%(mass/mass)
Density (20°C)	0.884~0.889g/ml
Refractive index 20°C	1.406~1.409
Water	max.0.05%(mass/mass)
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max.0.005%(mass/mass)
Cadmium (Cd) (RoHS)	max.10ppm
Mercury (Hg) (RoHS)	max.10ppm
Lead (Pb) (RoHS)	max.100ppm
Chromium (Cr) (RoHS)	max.10ppm
Brominated flame retardants (as Br) (RoHS)	max.1ppm