

Code No.20-1716

# SPECIFICATION

## L(-)-Threonine Standard

for Amino Acid Analysis

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
Specific rotation [ $\alpha$ ]D20(c=6,H2O)	-27.5~-29.0°
Loss on drying at 105°C	max.0.2%
Suitability for amino acids determination	under investigation
Assay (volumetric analysis)	min.99.0%
Assay (subtracting method)	under investigation