

Code No.20-1889

SPECIFICATION

L(+)-Tartaric Acid(Synthetic)

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
Specific rotation [α]D20(c=20,H2O)	+11.5~+13.5°
Residue after ignition (as sulfate)	max.0.1%
Chloride (Cl)	max.0.002%
Heavy metals (as Pb)	max.0.001%
Iron (Fe)	max.0.001%
Assay	min.99.0%