

SPECIFICATION

Fluorescein Sodium

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Solubility in water Ethanol-insoluble matter Absorbance (3mg/L, pH7.0) Loss on drying at 105°C Residue after ignition (as sulfate) (after drying) Chloride (Cl) Zinc (Zn) Sensitivity	Green or yellowish red ~ reddish brown, crystalline powder ~ powder or mass to pass test max.1.0%(mass/mass) min.0.45 max.7.0%(mass/mass) 35.0~40.5%(mass/mass) max.0.1%(mass/mass) max.0.05%(mass/mass) to pass test

JIS K 8830 : 2023(*:Additional test performed by Wako but not required by JIS)