

Code No.22-0088

## SPECIFICATION

## Vinyl Benzoate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ nearly colorless, clear liquid
Solubility in ethanol	to pass test
Density (20°C)	1.067~1.074g/ml
Refractive index 20°C	1.526~1.530
Water	max.0.3%
Acidity (as C <sub>6</sub> H <sub>5</sub> COOH)	max.0.2%
Assay (capillary GC)	min.98.0%