

Code No.22-0203

SPECIFICATION

Vascular Endothelial Growth Factor-A165, Mouse, recombinant, Animal-derived-free
for Cellbiology

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Bacterial endotoxins	less than 0.01 ng/ μ g
Electrophoresis test (SDS-PAGE)	to pass test
Mass spectrometry (MALDI-TOF)	to pass test

Source:E. coli expressed mouse vascular endothelial growth factor-A165
This product is recombinant cytokine which is produced in the process of culture and purification without using any substances from animals.