

## SPECIFICATION

Wakogel(R) C-200  
for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White small granules
Particle size (75~150 $\mu$ m)	min.75%
Particle size (pass 45 $\mu$ m)	max.5%
Precipitation volume	1.8~2.7ml/g
Loss on drying at 180°C	max.5.0%
pH (50g/l aqueous slurry,25°C)	5.5~7.0
Flowing rate	1.6~2.8ml/min
Separation power	to pass test
Activity	to pass test