

SPECIFICATION

Wakogel(R) C-300
for Column Chromatography

| REQUIREMENT | SPECIFICATION |
|---|---------------|
| Appearance | White powder |
| Particle size (45~75 μ m) | min.75% |
| Particle size (coarser than 75 μ m) | max.5% |
| Precipitation volume | 2.0~3.0ml/g |
| Loss on drying at 180°C | max.5.0% |
| pH (50g/l aqueous slurry,25°C) | 5.5~7.5 |
| Flowing rate | 0.5~1.4ml/min |
| Separation power | to pass test |
| Activity | to pass test |