

Code No.23-0143

SPECIFICATION

Wakogel(R) C-200E
for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White,
Particle size (75~150 μ m)	crystalline powder or small granules
Particle size (pass 45 μ m)	min.75%
Precipitation volume	max.5%
Loss on drying at 180°C	1.8~2.7ml/g
pH (50g/l aqueous slurry,25°C)	max.6.0%
Flowing rate	5.5~7.5
Separation power	1.6~3.0ml/min
Activity	to pass test
	to pass test