

SPECIFICATION

WakosilR C-200

for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White small granules
Particle size (63~212 μ m)	min.80%
Particle size (coarser than 212 μ m)	max.5%
Particle size (finer than 63 μ m)	max.10%
Precipitation volume	1.7~2.3mL/g
Loss on drying at 180°C	5.0~10.0%
pH(50g/L aqueous slurry, 25°C)	6.5~7.5
Flowing rate	1.6~2.3mL/min
Separation power	to pass test
Activity	to pass test