

Code No.23-0245

# SPECIFICATION

## WakosilR HC-N for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White powder
Particle shape	Spherical
Particle size (35~63 $\mu\text{m}$ )	under investigation
Particle size (coarser than 63 $\mu\text{m}$ )	under investigation
Particle size (finer than 35 $\mu\text{m}$ )	under investigation
Loss on drying at 180°C	under investigation
pH (50g/L aqueous slurry)	under investigation
Flow test	to pass test
Practical test	to pass test