

Code No.23-0270

## SPECIFICATION

WakogelR 60N, 63~212  $\mu$ m  
for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White small granules
Particle size (63~212 $\mu$ m)	min.75%
Particle size (pass 38 $\mu$ m)	max.10%
Precipitation volume	2.1~3.0mL/g
Solubility in ethyl acetate	to pass test
Loss on drying at 180°C	max.7.0%
pH(50g/L aqueous slurry, 25°C)	6.0~7.5
Flowing rate	1.3~3.2mL/min
Separation power	to pass test
Activity	to pass test