

SPECIFICATION

Xylene

Extra Pure

REQUIREMENT	SPECIFICATION
* Description Assay (o-, m-, p-Xylene) (GC) Appearance Density (20°C) Water Residue after evaporation Acidity (as HCl) Alkalinity (as NaOH) Sulfur compounds Thiophene group Benzene (GC) Toluene (GC) Substances darkened by sulfuric acid	Colorless clear liquid min.80%(mass/mass) max.10 Hazen units 0.860~0.870g/ml max.0.03%(mass/mass) max.0.002%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) to pass test to pass test max.0.1%(mass/mass) max.1.0%(mass/mass) to pass test

JIS K 8271:2007(*:Additional test performed by Wako)