

Code No.26-0013

## SPECIFICATION

## Zinc

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Bluish white or silver gray, mossy
Copper (Cu)	max.0.005%(mass/mass)
Cadmium (Cd)	max.0.005%(mass/mass)
Lead (Pb)	max.0.007%(mass/mass)
Arsenic (As)	max.0.1 ppm(mass/mass)
Iron (Fe)	max.0.005%(mass/mass)

JIS K 8012:2006

\*:Additional test performed by Wako but not required by JIS

This product is mossy.