

Code No.26-0013

## SPECIFICATION

## Zinc

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Copper (Cu) Cadmium (Cd) Lead (Pb) Arsenic (As) Iron (Fe)	Bluish white or silver gray, mossy max.0.005%(mass/mass) max.0.005%(mass/mass) max.0.007%(mass/mass) max.0.1 ppm(mass/mass) max.0.005%(mass/mass)

JIS K 8012:2006

\*:Additional test performed by Wako but not required by JIS  
This product is mossy.