

SPECIFICATION

Zinc

Reference Material for Volumetric Analysis

REQUIREMENT	SPECIFICATION
* Appearance	Bluish white turnings
Assay	min.99.99%(mass/mass)
Copper (Cu)	max.0.0010%(mass/mass)
Cadmium (Cd)	max.0.0010%(mass/mass)
Tin (Sn)	max.5ppm(mass/mass)
Lead (Pb)	max.0.0020%(mass/mass)
Iron (Fe)	max.0.0010%(mass/mass)

JIS K 8005:2006

*:Additional test performed by Wako but not required by JIS