

Code No.29-4745

SPECIFICATION

PresepR-C C18(ODS) (Short)

for Sample Pretreatment

| REQUIREMENT | SPECIFICATION |
|---------------------------|---------------------|
| Appearance | to pass test |
| Loss on drying at 85°C | max.5.0% |
| Weight of column packing | under investigation |
| Permitted pressure (1) | to pass test |
| Permitted pressure (2) | to pass test |
| Leakage of column packing | to pass test |
| Recovery | 90~105% |
| State of column packing | to pass test |