

SPECIFICATION

Diammonium Hydrogenphosphate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay [(NH ₄) ₂ HPO ₄] Solubility in water pH of a 50g/L solution at 25°C Chloride (Cl) Nitrate (NO ₃) Sulfate (SO ₄) Sodium (Na) Potassium (K) Arsenic (As) Lead (Pb) Iron (Fe)	White, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test 7.8~8.3 max.5ppm(mass/mass) max.0.003%(mass/mass) max.0.005%(mass/mass) max.0.01%(mass/mass) max.0.005%(mass/mass) max.1ppm(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass)

JIS K 9016:2023(*:Additional test performed by Wako but not required by JIS)