

SPECIFICATION

Potassium Hydrogen (+)-Tartrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (C ₄ H ₅ KO ₆) Solubility in dil.ammonia solution Chloride (Cl) Phosphate (PO ₄) Sulfur compounds (as SO ₄) Sodium (Na) Copper (Cu) Calcium (Ca) Zinc (Zn) Lead (Pb) Iron (Fe) Ammonium (NH ₄)	White, crystalline powder ~ powder min.99.5%(mass/mass) to pass test max.0.002%(mass/mass) max.0.002%(mass/mass) max.0.005%(mass/mass) max.0.1%(mass/mass) max.1ppm(mass/mass) max.0.005%(mass/mass) max.2ppm(mass/mass) max.1ppm(mass/mass) max.5ppm(mass/mass) max.0.002%(mass/mass)

JIS K 8537:2015

(*:Additional test performed by Wako but not required by JIS)