

SPECIFICATION

Potassium Dihydrogenphosphate for Oxidant Analysis
for Oxidant Analysis

REQUIREMENT	SPECIFICATION
* Appearance	White crystals
Assay (KH ₂ PO ₄)	min.99.0%(mass/mass)
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	4.2~4.6
Nitrogen compounds (as N)	max.0.001%(mass/mass)
Lead (Pb)	max.0.001%(mass/mass)
Suitability for oxidant determination	to pass test

JIS K 8118:2021(*:Additional test performed by Wako but not required by JIS)