

Code No.16-0939

## SPECIFICATION

### Potassium Dihydrogenphosphate for Oxidant Analysis for Oxidant Analysis

REQUIREMENT	SPECIFICATION
* Appearance	White crystals
Assay (KH <sub>2</sub> PO <sub>4</sub> )	min.99.0%(mass/mass)
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	4.2~4.6
Nitrogen compounds (as N)	max.0.001%(mass/mass)
Lead (Pb)	max.0.001%(mass/mass)
Suitability for oxidant determination	to pass test

JIS K 8118:2021(\*:Additional test performed by Wako but not required by JIS)