

Code No.18-0341

SPECIFICATION

Rutin Standard  
for Food Additives Test

REQUIREMENT	SPECIFICATION
Description	Pale yellow ~ pale yellowish green,crystalline powder
Identification	to pass test
Specific absorbance E1%1cm (λ max about 350nm)	min.290
Purity	-
Related substances	to pass test
* Solubility in methanol	to pass test
* Water	6.9~9.1%

This product is suited for Rutin for assay of reagents and test solutions on Japan's Specifications and Standards for Food Additives(\*:Additional test performed by Wako)