

Code No.18-0341

# SPECIFICATION

## Rutin Standard

### for Food Additives Test

REQUIREMENT	SPECIFICATION	
Description	Pale yellow ~ pale yellowish green, crystalline powder	
Identification		to pass test
Specific absorbance E1%1cm ( $\lambda$ max about 350nm)		min.290
Purity	-	
Related substances	to pass test	
* Solubility in methanol	to pass test	
* Water	6.9~9.1%	

This product is suited for Rutin for assay of reagents and test solutions on Japan's Specifications and Standards for Food Additives(\*:Additional test performed by Wako)