

SPECIFICATION

Tin(II) Chloride Dihydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (SnCl ₂ ·2H ₂ O) Solubility in hydrochloric acid Sulfate (SO ₄) Sodium (Na) Potassium (K) Copper (Cu) Calcium (Ca) Lead (Pb) Arsenic (As) Iron (Fe) Ammonium (NH ₄)	White crystals min.97.0%(mass/mass) to pass test max.0.002%(mass/mass) max.0.005%(mass/mass) max.0.002%(mass/mass) max.0.001%(mass/mass) max.0.002%(mass/mass) max.0.003%(mass/mass) max.2ppm(mass/mass) max.0.002%(mass/mass) max.0.002%(mass/mass)

JIS K 8136:2017*:Additional test performed by Wako but not required by JIS